

Veatch Software Specifications

	ReSeeVit Evolution	ReSeeVit 6.0	Huvitz HDRMate	Huvitz HIS-5000	HDC-7000
Processor	Pentium IV, 2 GHz - RAM: 2 GiG	Pentium IV, 2 GHz - RAM: 2 GiG	Equal or higher than 486 CPU with higher clock frequency than 90MHz, Main Memory: Minimum 8MB., HD Space: Minimum 2MB	CPU Pentium IV, 2GHz	CPU Pentium IV, 2GHz
Video card	ATI Radeon 9200 (128 MB) or equivalent	ATI Radeon 9200 (128 MB) or equivalent		Video Card ATI Radeon 9200 (128MB) or similar	
Computer interface	Firewire DCAM v1.30 400 Mb/s (6 Pin) - No Built In Motherboard	USB	USB	Firewire DCAM v1.30 400 Mb/s (6 Pin) - No Built In Motherboard	USB
Color monitor	LCD or CRT, 1280 x 1024 resolution @75 Hz (CRT with 0.24 mm pitch)	LCD or CRT, 1280 x 1024 resolution @75 Hz (CRT with 0.24 mm pitch)			LCD Compatibility - 17" (1280x1024) 19"(1280x1024) 20.1"(1400x1050) 20.1"(1680x1050) 21"(1680x1050) 22"(1680x1050) 20.1"(1600x1200) 23"(1920x1080) 23"(1920x1200) 24"(1920x1200) 25"(1920x1080)
Operating systems:	Windows 7 Pro (32 Bit), or Windows XP Pro operating system. Note: If Windows XP is installed, you will need to install Service Pack 2. Internet Explorer 6 must be installed.	Windows 7 Pro (32 Bit), or Windows XP Pro operating system. Note: If Windows XP is installed, you will need to install Service Pack 2. Internet Explorer 6 must be installed.	Windows 95/98/ME, 2000, XP,7 (32 & 64 Bit) OTHER Software Components: Microsoft Office 97 or higher.	Windows 7 Pro (32 & 64 Bit), or Windows XP Pro operating system. Note: If Windows XP is installed, you will need to install Service Pack 2. Internet Explorer 6 must be installed.	Windows 2000, XP, Vista, 7 (32 & 64 Bit)

Sirius Acuity Software

E series Intel Pentium or equivalent Athlon processor (Performance will be best on any iCore Intel eg: i3, i5, or i7);Hard Drive 300+ GB 2 GB RAM
Dedicated graphics card: ATI Radeon Series or nVidia GeForce Series
USB
Minimum Size: 17" Diagonal Monitor
Windows XP or Windows 7, (32 bit)