

STANDARD KERATOMETER

RIGHTmed
S4OPTIK



S4OPTIK


Obtain accurate corneal measurements quickly and easily with the S4OPTIK Keratometer. Built-in contact lens holder also allows measurement of rigid contact lenses.

TECHNOLOGY
VALUE
RELIABILITY

FEATURES:

- ◆ coincidence focusing system. Positive measurements to 0.02mm accuracy
- ◆ one-position instrument. Measure both meridians without changing optical system
- ◆ allows qualitative review of corneal surface
- ◆ measurement dials provide both diopter and millimeter values
- ◆ sturdy construction for extended life
- ◆ fits all standard instrument trays

Clear mires make alignment and focus effortless



Fast, Easy Corneal Curvature Measurements



SPECIFICATIONS

Measuring range	Radius of curvature, 9.4 ~ 6.4 mm (In 0.05 mm steps), Diopter 36 ~ 52 (in 0.25 steps).
Instrument Dimensions	Height - 16", Width - 7", Length - 19"
Dimensions	16" H (406 mm) x 19" L (483 mm) x 7" W (178 mm)
Weight of instrument	7.0 Kgs approximately
Colour / Finish	Black
Mode of operation	Continuous
Lamp	230v-15w, 110v-25w, E12 Type
Standard Accessories	Dust cover, Spare Bulbs (2 No's)
Optional Accessories	Stilts for table mount Tray for stand mount

Javal Ophthalmometer also available from S4OPTIK



Distributed by:

with a policy of continuous improvement, we reserve the right to make changes to specifications without advance notice.

RIGHTmed
S4OPTIK

5325 Cleveland Street, Suite 303 Virginia Beach, VA 23462
 1.888.224.6012 | www.s4optik.com